



 **Quaha Zono**  
Enabling | Independent | Life



# QUHA ZONO IS THE ONLY WIRELESS GYROSCOPIC MOUSE DESIGNED FOR SPECIAL NEEDS.

It enables independent computer access.

## IS QUHA ZONO FOR ME?

> It provides accurate and intuitive computer access with just small head movements!

## CAN I USE QUHA ZONO WITH MY DEVICES?

> It is compatible with all available computers and almost all tablets and smartphones! \*)

## HOW DO I START USING QUHA ZONO?

> There is absolutely zero setup time: open the box, plug it in and turn it on. Instant access to your device!

\*) currently limited compatibility with Apple iOS devices (iPhone & iPad)



Quha Zono mouse  
USB receiver  
Universal clip  
USB charging cable  
Adapter to connect two switches

User manual and CD  
Warranty card  
Zono button adapter  
Quha velcro band

## TECHNICAL SPECIFICATIONS

|                     |                                                               |
|---------------------|---------------------------------------------------------------|
| Dimensions          | Length 59 mm, width 33 mm, height 15 mm and weight 23 g       |
| Wireless connection | Radiolink up to 10 m                                          |
| Battery             | Active use time 30 h, in automatic idle mode up to 20 days    |
| USB receiver        | USB 2.0 HID                                                   |
| Switch connectivity | 3.5 mm connector on both Quha Zono mouse and the USB receiver |

# Quha Zono enables you to be you.

Quha Zono makes it possible to use a computer in a flexible way, when you cannot use an ordinary mouse. It is an extremely light device, and since it is completely wireless there are endless ways to wear it.

Quha Zono works straight out of the box allowing full access to your computer. Take complete control with wide range of accessories. Quha Zono represents functional Nordic design and is designed and manufactured in Finland.



I have been waiting for this mouse for 10 years.

- Paula, Zono user from Finland



Designed and manufactured in Finland by Quha oy

[www.quha.com](http://www.quha.com)

Copyright 2014 Quha oy.

Quha, the Q-mark and Quha Zono are trademarks or registered trademarks of Quha oy