

ASSISTIVE TECHNOLOGY PRODUCT CATALOG 2020



QUHA ZONO - GYROSCOPIC WIRELESS MOUSE



Quha Zono is the original gyrosopic mouse designed for special needs.

Quha Zono makes it possible to use a computer when you cannot use an ordinary mouse with your hands. It translates natural head movements into mouse pointer movements – the mouse pointer follows your head movements.

It is an extremely light device, and since it is completely wireless there are endless ways to wear it.

WITH QUHA ZONO HEAD MOUSE YOU CAN DO EVERYTHING YOU CAN DO WITH THE MOST ACCURATE DESKTOP MICE.

Quha Zono is compatible with all Windows, macOS, Chromebook and Linux computers, iPads, Android tablets and most smartphones.

Sales package includes:

- Quha Zono mouse
- USB Receiver
- USB charging cable
- Stereo button adapter for Zono
- Adapter to connect two switches
- CD with configuration application
- User manual and warranty card

Choose the best way to wear Quha Zono from the options presented on page 4.



//

I have been waiting for this mouse for 10 years.

P/N	PRODUCT
20	QUHA ZONO GYROSPIC MOUSE

QUHA ZONO 2 - ENHANCED GYROSCOPIC WIRELESS MOUSE

Quha Zono 2 is an enhanced version of the legendary Quha Zono gyroscopic mouse. It includes several new features requested by the users, for example, cursor centering, pausing and scrolling. These can be activated either by simple gestures or with buttons connected to Zono or USB Receiver.

In addition, the active usage time of Quha Zono 2 is over doubled to 60 hours! It also has support for smart accessories.

QUHA ZONO 2 INCLUDES INTEGRATED CURSOR CENTERING, PAUSING AND SCROLLING FEATURES.

Sales package includes:

- Quha Zono mouse
- USB Receiver
- [HEAD MOUNT KIT \(SEE PAGE 4\)](#)
- [LICENSE FOR QUHA DWELL \(PAGE 9\)](#)
- USB charging cable
- Stereo button adapter for Zono
- Adapter to connect two switches
- CD with configuration application
- User manual and warranty card

Quha Zono 2 is compatible with all current Zono accessories. Choose the best way to wear Quha Zono 2 from the options presented on page 4.



//
Love to have a single receiver with
which you can choose if the
switch is click, center or scroll!

P/N	PRODUCT
21	QUHA ZONO 2 GYROSCOPIC MOUSE

HEADBAND

Headband is an easy way to wear Quha Zono. It is light and fits comfortably for hours. Headband is adjustable to make it suitable for most users.

Headband is an ideal alternative when using Zono while resting the head on a pillow or a headrest.

Quha Zono is easily attached directly to Headband, making it very light to wear. The Quha Zono mouse can be charged while attached to Headband.



EYEWEAR KIT

Eyewear kit is a spectacle type way to wear Quha Zono. It includes a light and firm eyeglass frame and the Light Eyeglass Clip.



UNIVERSAL WEAR KIT

Universal Wear Kit enables Quha Zono usage by hand or by foot. It includes a Velcro band and a soft holder for Quha Zono.



HEAD MOUNT KIT

Head Mount Kit fits comfortably on most users' heads and can be worn for hours.



BASEBALL CAP

Quha Baseball cap includes a classic cap and a white soft holder for Quha Zono.



LIGHT EYEGLASS CLIP

Light Eyeglass clip is an ultra light clip to attach Quha Zono safely to your own eyeglasses.



P/N	PRODUCT
170	HEADBAND, ADJUSTABLE
101	HEAD MOUNT KIT
103	LIGHT EYEGLASS CLIP
130	EYEWEAR KIT
111	UNIVERSAL WEAR KIT
132	BASEBALL CAP

QUHA ZONO USER STORIES

“IT IS VERY IMPORTANT THAT I CAN USE A COMPUTER AND A CAMERA BY MYSELF”

My name is Jari Mönkkönen and I am tetraplegic. I worked as a firefighter for 12 years. I also have a long career in artistic gymnastics. I have 29 gold medals at the Finnish championships and also a silver medal at the European championship on high-bar in '98.

My tetraplegia is caused by a sports accident. I made a double backflip and it landed on my neck. Nowadays I work as a head coach in a company called Voltimedia. We train young people and sports clubs.

Computers are very important in my daily life. After this accident, in a rehabilitation center, I tried several different devices and chose the Quha Zono head mouse. Quha Zono is accurate, fast, light-weight and easy to use.



“I WOULD DEFINITELY RECOMMEND THIS DEVICE TO OTHERS”



My name is Katherine Quinn. I am 18 years old and I have cerebral palsy. My muscles don't work right, but my brain works just fine. I live in Carmichael, California. My interests are scrapbooking, horseback riding and watching videos. I'm going to be a senior at Rio Americano High school, where my favorite class is music. One day I would like to be Judge Katherine Quinn.

I have previously tried a headmouse, but it was difficult to use. I also tried the Eye Max, but it was too hard to set up every time I changed classes. Then I heard about the Quha products from Maggie Mahoney at Quha USA.

With the Quha Zono mouse I can easily upload a video on YouTube and read ebooks all by myself. Another great thing about the Quha Zono mouse is that you can move around with it – this is very important for me as sometimes my body moves by itself and I can't control it. These are the reasons why Zono is so great. I would definitely recommend this device to others.

See full user stories at www.quha.com/stories

QUHA NEMO - SELF OPERATED HEADBAND WITH CHARGING DOCK



Quha Nemo is a modified version of the Headband. It includes a charging dock with a possibility to pick it up independently without any assistance.

Quha Nemo can also be used as an independent charging dock with easy attachment without disconnecting cables.

Nemo is designed so that the user can pick up the Headband with their head and start using it immediately. There is no need for assistance, and after using it, the user can return it to the charging dock where it will automatically start charging.

There are two variants of the Nemo, one with a Sento puff switch and one without. Quha Nemo comes with a Zono mouse, Nemo Headband, optional Sento Puff switch, the charging dock and a charger.



INDEPENDENT WAY TO WEAR AND USE THE ZONO MOUSE.

WHO IS IT FOR?

Nemo is for those who wish to use the Zono mouse independently without any assistance.

HOW DO I INSTALL THE NEMO?

Installation is very easy. All you need is a mount with a camera thread. Attach the Nemo to it, plug in the charger and you're good to go.

P/N	PRODUCT
80	QUHA ZONO WITH NEMO
81	QUHA ZONO WITH NEMO AND SENTO

Quha Nemo makes it possible to use Quha Zono independently.

QUHA SENTO - CONTACTLESS PUFF SWITCH



Quha Sento is the world's first contactless puff switch. It is used to perform a short mouse click or dragging together with the Quha Zono gyroscopic mouse.

Mouse clicks are made easily and hygienically with a small puff to the microphone type of mouthpiece.

WITH QUHA SENTO YOU CAN PERFORM MOUSE CLICKS EASILY, HYGIENICALLY AND WITHOUT HANDS.

DESIGNED CONSUMER ELECTRONICS LOOK & FEEL

- No internal mouth pieces

COMPLETE WIRELESS COMPUTER ACCESS


- Sento together with Quha Zono is the complete solution for computer access.

THE ONLY 100% HYGIENIC PUFF SWITCH ON THE MARKET

- Enables one device to be shared between multiple users.

NO CONTROL DELAY, UNIQUE SENSOR TECHNOLOGY

- Instant puff response



//

The design and the feeling of quality are great!

P/N	PRODUCT
190	QUHA SENTO WITH ADJUSTABLE HEADBAND
191	QUHA SENTO WITH HEAD MOUNT KIT
192	QUHA SENTO WITH EYEWEAR KIT

QUHA PUFO+ - PUFF SWITCH FOR QUHA ZONO

Quha Pufo+ is a wireless puff switch for the Quha Zono gyroscope mouse, making it even more versatile as a head mouse.

Quha Pufo+ puff switch is quickly connected to Quha Zono — no additional power or configurations are needed. The mouthpiece is easily adjustable to each user and replaceable for hygiene.

A little puff to Pufo's mouthpiece will generate a left mouse click and two short puffs a double-click. The dragging can be easily



P/N	PRODUCT
181	QUHA PUFO+ WITH HEADBAND
182	QUHA PUFO+ WITH HEAD MOUNT KIT

QUHA SENTO MP - MULTIPURPOSE CONTACTLESS PUFF SWITCH

Quha Sento MP is a generic multipurpose contactless puff switch. Quha Sento MP enables users to control any assistive technology system in the world. For example, environmental control, eye tracking, nurse calls and other systems can easily be controlled by Quha Sento MP. It replaces any ability switch with the same standard 3.5 mm connector.

Quha Sento MP includes a flexible sensor arm, a camera thread, a buzzer and configurable sensitivity.



P/N	PRODUCT
193	QUHA SENTO MP

QUHA DWELL - DWELL CONTROL FOR COMPLETE COMPUTER ACCESS

Quha Dwell is a complete computer access software for assistive mice.

The unique Instant Dwell feature makes computer usage faster, more comfortable and ergonomic compared to any current dwell control software.

With Quha Dwell you can easily control all the functions of your computer, switch between applications, copy-paste and scroll. Any new function can be added to Quha Dwell to meet your needs.

Quha Dwell includes an on-screen-keyboard with word prediction.



P/N	PRODUCT
270	QUHA DWELL

OPTIONS TO CONNECT QUHA ZONO TO AN IPAD OR IPHONE

Quha Zono is connected to an iPad or iPhone with a Lightning adapter and to an iPad Pro with a USB-C adapter. In both cases Quha Zono is a plug-and-play device and works without any additional software.

Users with an older iPad/iPhone that cannot be updated to iPadOS/iOS 13 can use an AMANEO BTi adapter.

P/N	PRODUCT
211	IPAD ADAPTER FOR LIGHTNING PORT
221	USB-C OTG ADAPTER
222	USB-C HUB
210	AMANEO BTi IPAD ADAPTER



QUHA KAJO BUTTON - ABILITY SWITCH



Quha Kajo Button is a sensitive and durable Assistive Technology switch designed for everyday use. It can be used in any situation that requires switch access.

Quha Kajo Button has been designed to endure hard use and no compromises have been made with the quality of components. It is mountable using the thread inserts in the bottom.

Quha Kajo Button switches have distinctive auditory and tactile feedback. They are available in two different sizes - 25 mm (1") and 50 mm (2").

QUHA KAJO SOFT BUTTON - PILLOW SWITCH



Quha Kajo Soft Button is a soft switch with a large soft surface. Because of the soft foam layer under the cover the Soft Button can be safely pressed with head movement.

The safety pin in the bottom of the Soft Button ensures that the button stays where it is supposed to stay during the night. It can be attached to a pillow or other fabrics. The safety pin is detachable.

The cover of the Soft Button is dirt resistant and skin friendly. The cover can be washed (up to 60°C/140°F)

P/N	PRODUCT
250	QUHA KAJO BUTTON 25 MM
251	QUHA KAJO BUTTON 50 MM
252	QUHA KAJO SOFT BUTTON

MORE INFORMATION
ON ALL QUHA
PRODUCTS:
www.quha.com/products



QUHA VIDO - EXTERNAL VIRTUAL KEYBOARD AND MOUSE

Quha Vido is an external virtual keyboard and mouse. Quha Vido allows you to use your computer more efficiently, since it frees up your screen space as there is no need to display a virtual keyboard and mouse on your computer screen. Using Quha Vido makes writing emails, browsing the internet and working with your computer faster and more efficient.

The easy and fast setup lets you decide where you want to place the Vido display. You can use Quha Vido together with any wireless computer mouse provided by other manufacturers, although using it in conjunction with the Quha Zono wireless computer mouse gives you the best user experience. With the help of its adjustable stand, you can place the Vido display exactly how you want it and adjust the angle to your liking.

EFFORTLESS DISPLAY FOR UNLIMITED COMPUTER ACCESS

- Plug-and-play.
- Automatic one-click installation for virtual keyboard and mouse.
- Place the screen where it best suits you.

MAKES IT POSSIBLE TO USE COMPUTER EVEN FASTER

- Make room for important things on screen by storing keyboard on an external screen.
- External USB monitor for virtual keyboard and mouse.
- Powered from a USB cable.

//

I love the Quha Vido. I can use my computer much faster now

Vido user from Finland



P/N

PRODUCT

300

QUHA VIDO



Designed in Finland by Quha

www.quha.com

Copyright 2020 Quha oy.

Quha, the Q-mark, Quha Zono, Quha Zono 2, Quha Sento, Quha Sento MP, Quha Pufo+, Quha Vido, Quha Dwell and Instant Dwell are trademarks or registered trademarks of Quha oy.

iPhone, iPad and AssistiveTouch are trademarks of Apple Inc., registered in the U.S. and other countries. iOS is a trademark or registered trademark of Cisco in the U.S. and other countries and is used under license.